

Whether sample exactly matches with the target?



The difference is small enough to be detected or accepted by





SHAKTI VINAY SHUKLA

FRAGARNCE & FLAVOUR DEVELOPMENT CENTRE, KANNAUJ



OLFACTION

Material:Smelling Strips/SlipsSize:Length = 14 - 18 cmWidth = 5-10 mmWeight:180 gsmfree from any odour



ATMOSPHERE FOR OLFACTION

- *** ODOURLESS**
- * HUMID
- * CALM
- * TEMPERATE
- *** FREE FROM DISTRACTION**



TIME & ORDER OF OLFACTION

- MORNING IN FRESH AIR
- FROM WEAK TO STRONG ODORANT
- NO RANDOM SMELLING
- SYSTEMATIC



REGENERATION & BREATHING

- BREAK BETWEEN EACH SERIES
- OXYGENATION BY BREATHING LEADS TO REGENERATION
- **BREATH IN THROUGH NOSE**
- BREATH OUT THROUGH MOUTH





- Marked one end before use
- Dip 0.5-1.0 cm in fragrant material
- No touch to nose
- **3-4** breath enough
- I-10% Solution for strong odorant



Test Form

Analytical

- 1. Discriminative
- 2. Descriptive



Affective

- **1. Paired Preference**
- 2. Hedonic rating
- 3. Ranking

Analytical Test Form

Discriminative

- **1.** Paired Comparison
- 2. Duo-trio test
- **3.** Triangle Test
- 4. Two out five test
- 5. A or not A test

Descriptive

- Simple
- Rating
- Odour Profile
- Dilution profile technique



Profile

